# **Impella CP® with SmartAssist®**

**USER MANUAL  
GEBRAUCHSANWEISUNG  
MANUEL D'UTILISATION  
MANUAL DEL USUARIO  
MANUALE UTENTE  
BRUKSANVISNING  
BRUGERVEJLEDNING  
GEBRUIKERSHANDLEIDING  
BRUKSANVISNING  
INSTRUKCJA OBSŁUGI  
UŽIVATELSKÁ PŘÍRUČKA  
MANUAL DO USUÁRIO**

**OO  
OO ABIOMED®**

Recovering hearts. Saving lives.

## TABLE OF CONTENTS

English .....	1
Deutsch .....	15
Français .....	29
Español.....	43
Italiano .....	57
Svenska.....	71
Dansk .....	85
Nederlands .....	99
Polski.....	113
Norsk.....	127
Česky.....	139
Portugues .....	151

## USER RESPONSIBILITY

The Impella System will perform in conformity with the description thereof contained in this manual and accompanying labels and/or inserts, when assembled, operated, and maintained in accordance with the instructions provided. A defective system should not be used. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## INTENDED AUDIENCE

This manual is intended for use by cardiologists, surgeons, nurses, perfusionists, and cardiac catheterization lab technicians who have been trained in the use of the Impella Circulatory Support System.

## WARNINGS

Warnings alert you to situations that can cause death or serious injury. The dark gray  symbol appears before warning messages.

	Use of the Impella System by trained and experienced practitioners has been associated with improved outcomes. Consequently, the first use of Impella should be preceded by the completion of a contemporary Abiomed Impella training program and include on-site proctoring during the first use by Abiomed clinical support personnel certified in the use of Impella.
	Fluoroscopy is required to guide placement of the Impella Catheter. The 0.018 inch placement guidewire must be reliably observed at all times.
	Be sure that the stopcock on the introducer or repositioning sheath is always kept in the closed position. Significant bleed back can result if the stopcock is open.
	Avoid manual compression of the inlet and outlet areas of the cannula assembly.
	The sterile components of the Impella System can be used only if the sterilization indicators show that the contents have been sterilized, the packaging is not damaged, and the expiration date has not elapsed.
	Do <b>NOT</b> resterilize or reuse the Impella Catheter. It is a disposable device and is intended for single use only. Reuse, reprocessing, or resterilization may compromise the structural integrity of the catheter and/or lead to catheter failure which, in turn, may result in patient injury, illness, or death.
	Retrograde flow will occur across the aortic valve if the Impella Catheter is set at P-level P-0.
	To prevent malfunction of the locking mechanism of the introducer, do <b>NOT</b> hold the hemostatic valve while inserting into the artery.
	To prevent failure of the peel-away introducer, remove the peel-away introducer prior to transport when activated clotting time (ACT) is 160-180 seconds.
	Do <b>NOT</b> use saline in the purge system.
	Do <b>NOT</b> use an Impella System if any part of the system is damaged.
	To prevent the risk of explosion, do <b>NOT</b> operate the Impella System near flammable anesthetics.
	If at any time during the course of support with the Impella Catheter, the Automated Impella Controller alarms "Purge Pressure Low" or "Purge System Open," follow the instructions presented the Automated Impella Controller User Manual

   
Do **NOT** subject a patient which is supported by an Impella Catheter to magnetic resonance imaging (MRI). The strong magnetic energy produced by an MRI machine may cause the Impella System components to stop working, and result in injuries to the patient. An MRI may also damage the electronics of the Impella System.

 Cardiopulmonary support (CPR) should be initiated immediately per hospital protocol if indicated for any patient supported by the Impella Catheter. When initiating CPR, reduce the Impella Catheter flow rate. When cardiac function has been restored, return flow rate to the previous level and assess placement signals on the controller.

 During defibrillation, do **NOT** touch the Impella Catheter, cables, or Automated Impella Controller.

 Infusion through the sideport of the introducer can be done only after all air is removed from the introducer. If performed, the infusion should be done for flushing purposes only and **NOT** for delivering therapy or monitoring blood pressure.

 In patients with transcatheter aortic valves position the Impella system carefully to avoid interaction with the transcatheter aortic valve prosthesis. Unintentional interaction of the Impella motor housing with the TAVR device may result in destruction of the impeller blades. This can lead to systemic embolization, serious injury, or death. In this situation, avoid repositioning while the device is running; turn the device to P0 during repositioning or any movement that could bring the outlet windows into proximity to the valve stent structures. If there is low flow observed in a patient implanted with a transcatheter aortic valve prosthesis, consider damage of the impeller and replace the Impella as soon as possible.

 To reduce the risk of cardiac or vascular injury (including perforation) when manipulating the heart during cardiac surgery, evaluate the position of the pump using imaging guidance prior to manipulating the heart, and monitor position.

 To reduce the risk of cardiac or vascular injury (including ventricular perforation) when advancing or torquing the Impella, adjustments should be performed under imaging guidance.

 To reduce the possibility of fibers being drawn into the Impella, customers should avoid exposing the inlet and cannula section of the Impella Heart Pumps to any surfaces or fluid baths where the device can come into contact with loose or floating fibers.

 To reduce the risk of cardiac injury (including ventricular perforation), physicians should exercise special care when inserting the Impella Catheter in patients with complex anatomy. This includes patients with known or suspected: decreased ventricular cavity size, ventricular aneurysms, congenital heart disease, or compromised cardiac tissue quality in the settings of acute infarction with tissue necrosis.

 To reduce the risk of vascular injury, physicians should exercise caution when inserting the Impella Catheter in patients with complex peripheral vascular anatomy. This includes patients with known or suspected: unrepaired abdominal aortic aneurysm, significant descending thoracic aortic aneurysm, dissection of the ascending/transverse/descending aorta, chronic anatomical changes in the relationship of the aorta/aortic valve/ventricular alignment, significant mobile atherosomatous disease in the thoracic or abdominal aorta or peripheral vessels.

 Physicians should exercise special care when inserting the Impella Catheter during active Cardiopulmonary Resuscitation (CPR). In addition, active CPR maneuvers may change the position of the Impella Device, introducing the risk of cardiac or vascular injury (including ventricular perforation). Check that the pump is positioned correctly after CPR with echocardiography guidance.

## CAUTIONS

Cautions indicate situations in which equipment may malfunction, be damaged, or cease to operate. The light gray  symbol appears before caution messages.

-  Handle with care. The Impella Catheter can be damaged during removal from packaging, preparation, insertion, and removal. Do **NOT** bend, pull, or place excess pressure on the catheter or mechanical components at any time.
-  Patients with aortic stenosis or other abnormal aortic valve performance may be compromised by the use of the Impella Catheter. Patients with aortic valve disease should be observed for aortic insufficiency.
-  Partial circulatory support with Impella has been associated with more extensive use of rotational atherectomy. Extensive use of rotational atherectomy has been associated with a periprocedural increase in cardiac biomarkers indicative of myocardial injury. Rotational atherectomy, with or without the use of hemodynamic support, should be used in accordance with the manufacturer's instructions for use.
-  Use only original accessories and replacement parts supplied by Abiomed.
-  To prevent device failure, do **NOT** start the Impella Catheter until the guidewire has been removed.
-  Do **NOT** remove the Impella Catheter over the length of the guidewire.
-  When replacing the purge cassette, the replacement process must be completed within 90 seconds. The Impella Catheter may be damaged if replacement takes longer than 90 seconds.
-  Do **NOT** kink or clamp the Impella Catheter or the introducer.
-  Have a backup Automated Impella Controller, purge cassette and Impella Catheter available in the unlikely event of a device failure.
-  Do **NOT** alter the Impella Introducer kit in any way.
-  Aspiration and saline flushing of the Impella Introducer kit sheath, dilator, and valve should be performed to help minimize the potential for air embolism and clot formation.
-  Indwelling introducer sheaths should be internally supported by a catheter or dilator.
-  Dilators and catheters should be removed slowly from the sheath. Rapid removal may damage the valve, resulting in blood flow through the valve.
-  Never advance the guidewire or sheath when resistance is met. Determine the cause of resistance using fluoroscopy and take remedial action.
-  When injecting or aspirating through the sheath, use the sideport only.
-  Operation of the system without heparin in the purge solution has not been tested. In the event that a patient is intolerant to heparin, due to heparin-induced thrombocytopenia or bleeding, physicians should use their clinical judgment to assess the risks versus benefits of operating the Impella System without heparin. If it is in the best interest of the patient to operate the system without heparin, the glucose solution is still required, and physicians should consider systemic delivery of an alternative anticoagulant. Do **NOT** add any alternative anticoagulant (such as a direct thrombin inhibitor) to the purge fluid. The Impella Catheter has not been tested with any alternative anticoagulants in the purge solution.
-  Do not insert the repositioning sheath into the peel-away sheath during catheter insertion; slide the repositioning sheath towards the red pump handle prior to inserting the device.
-  Never remove the guidewire when increased resistance is met. Determine the cause of resistance using fluoroscopy and verify the repositioning sheath is not inserted into the peel-away sheath.
-  Use with Shockwave Coronary Intravascular Lithotripsy (IVL) Catheter at a distance of less than 20 mm between the optical sensor and IVL device may interfere with or damage the Impella optical sensor. Prior to IVL therapy pulsing, physicians should assess and verify this distance. If the placement signal is not displayed, monitor patient hemodynamics and confirm Impella position with imaging and motor current pulsatility. Loss of the placement signal does not impact Impella hemodynamic support.

## INTENDED PURPOSE

The Impella (intracardiac pump for supporting the left ventricle) is intended for clinical use in cardiology and in cardiac surgery for up to 5 days for the following indications, as well as others:

- The Impella is a circulatory support system for patients with reduced left ventricular function, eg, post-cardiotomy, low output syndrome, cardiogenic shock after acute myocardial infarction, or for myocardial protection after acute myocardial infarction
- The Impella may also be used as a cardiovascular support system during coronary bypass surgery on the beating heart, particularly in patients with limited preoperative ejection fraction with a high risk of postoperative low output syndrome
- Support during high risk percutaneous coronary intervention (PCI)
- Post PCI

## CONTRAINDICATIONS

- Mechanical aortic valves, severe aortic valvular stenosis or valvular regurgitation
- Hematological disorder causing fragility of the blood cells or hemolysis
- Hypertrophic obstructive cardiomyopathy (HOCM)
- Aneurysm or necrotomy or severe anomaly of the ascending aorta and/or the aortic arch
- Mural thrombus in the left ventricle
- Ventricular septal defect (VSD) after myocardial infarction
- Anatomic conditions precluding insertion of the pump
- Other illnesses or therapy requirements precluding use of the pump
- Severe peripheral arterial occlusion disease (PAOD) is a relative contraindication

## POSSIBLE COMPLICATIONS

There are risks of complications with every procedure using a blood pump. These include among others:

- Acute renal dysfunction
- Aortic valve injury
- Cardiogenic shock
- Cerebral vascular accident/Stroke
- Myocardial infarction
- Renal failure
- Thrombocytopenia
- Hemolysis
- Bleeding
- Limb ischemia
- Immune reaction
- Embolism, thrombosis
- Cardiac or vascular injury (including ventricular perforation)
- Positioning problems causing haemolysis or reduced haemodynamic support
- Infection and septicemia
- Dislocation of the pump
- Cardiovalvular injuries due to extreme movement of the suction cannula in relation to the cardiac valve or as a result of attachment by suction of the pump to the valve system following incorrect positioning
- Endocardial injuries as a result of attachment of the pump due to suction
- Pump failure, loss of pump components following a defect
- Patient dependency on the pump after use for support

## OVERVIEW

The Impella Catheter is an intravascular microaxial blood pump that supports a patient's circulatory system. The Impella Catheter is inserted percutaneously through the femoral artery and into the left ventricle.

When properly positioned, the Impella Catheter delivers blood from the inlet area, which sits inside the left ventricle, through the cannula, to the outlet opening in the ascending aorta. Physicians and device operators monitor the correct positioning and functioning of the Impella Catheter on the display screen of the Automated Impella Controller.

The clinical benefit of the Impella CP with SmartAssist is described in the Summary of Safety and Clinical Performance (SSCP) and can be accessed in EUDAMED (when operational) or requested from Abiomed Europe GmbH.

This section describes the components of the Impella Catheter and the Automated Impella Controller, as well as the accessory components.

### REUSABLE SYSTEM COMPONENTS

The Impella System consists of the following reusable components:

- Automated Impella Controller—provides the user interface, alarm indications, and portable battery
- Automated Impella Controller cart—for easy transport of the Automated Impella Controller

### SINGLE-USE SYSTEM COMPONENTS

The Impella System also includes the following single-use components:

- Impella Catheter
- Purge cassette
- Introducer kit
- 0.018 inch, 260 cm placement guidewire

### IMPELLA SET-UP AND INSERTION KIT

The Impella Catheter and the Single-Use components provided by Abiomed for the set-up and insertion of the Impella Catheter are provided in the Impella Set.

The Impella CP® with SmartAssist® Set contains:

- Impella CP® with SmartAssist® Catheter
- 0.018 inch, 260 cm placement guidewire
- Purge cassette
- Introducer kit
  - » Peel-away introducers - two 14 Fr introducers (13 cm and 25 cm)
  - » Dilators - 8 Fr, 10 Fr, 12 Fr, and 14 Fr
  - » 0.035 inch stiff access guidewire

### SYSTEM CONFIGURATIONS

The Automated Impella Controller connects to the Impella Catheter and accessory components in the initial set-up configuration as illustrated.


Figure 1.1

## IMPELLA® CATHETER

The Impella Catheter is an intravascular microaxial blood pump. The Impella CP® with SmartAssist delivers up to a maximum flow rate of 4.3 liters of blood per minute from the left ventricle into the aorta.


Figure 1.2

## PURGE CASSETTE

Do not use saline in the purge system.

The purge cassette delivers rinsing fluid to the Impella Catheter. The purge fluid (typically 5% glucose solution in water with heparin) flows from the purge cassette through the catheter to the microaxial blood pump to prevent blood from entering the motor. When the purge cassette is properly installed in the Automated Impella Controller, the Abiomed® logo is upright and facing you.

If included, remove and discard the Y Connector from the tubing when setting up the Impella CP® with SmartAssist®. Connect the yellow luer on the end of the purge tubing directly to the yellow luer on the Impella Catheter.


Figure 1.3

## ACCESSORIES

### Component Descriptions

The Impella CP® with SmartAssist® introducer kit is used to gain arterial access for the Impella CP® with SmartAssist® Catheter. It contains:

- Two 14 Fr (13 cm and 25 cm) peel-away introducers—with hemostatic valve for tight fit around components and single-step "break-away" configuration
- 8 Fr, 10 Fr, 12 Fr, and 14 Fr dilators—easy to insert and remove with soft design for atraumatic approach into femoral artery
- 0.035 inch stiff access guidewire

The 0.018 inch, 260 cm placement guidewire is used for the placement of the catheter. The guidewire has a radiopaque, shapable tip.

*These alternative guidewires have been tested and approved for use with the Impella System: Boston Scientific Platinum Plus™ ST 0.018 in, and Boston Scientific V-18 Control Wire™ ST 0.018 in.*

Hospital Provided:

Glucose solution (typically 5% glucose in water with 50 IU/mL of heparin) is used as the purge fluid through the Impella Catheter.

The Automated Impella Controller cart holds the Automated Impella Controller. The cart has wheels for easy transport of the controller and a storage basket. (For more information, including assembly instructions, refer to the Automated Impella Controller cart instructions for use.)

## PRE-SUPPORT EVALUATION

Before initiating the procedure, evaluate the patient for factors that may prevent successful placement of the Impella Catheter. Use imaging technology to examine the patient's vasculature and femoral access site. An echo assessment of the left ventricle is also recommended to rule out left ventricular thrombus, mechanical aortic valves, or severe aortic insufficiency. If the patient has a tortuous iliac artery, consider using the provided 25cm sheath for insertion. The Impella Catheter can also be inserted surgically.

## STARTUP

### SUPPLIES NEEDED

- Automated Impella Controller (Abiomed)
- Impella Set (Abiomed)
- Hospital provided:
  - Diagnostic catheter (AL1 or MP without side holes or pigtail with or without side holes)
  - Standard 0.035" x 175 cm J-tip guidewire
  - 500 cc bag of glucose solution for purge solution (5% recommended; 5% to 20% acceptable) with 50 IU heparin/mL

**Note:** If the purge solution is supplied in bottles, open the vent on the purge fluid spike and follow the same procedure as if supplied in bags.

## CASE START

To avoid fibers drawn into the Impella:

- Keep the Impella Heart Pump in its packaging tray until just before insertion.
- Do not attempt to run the pump in a basin of saline prior to insertion.
- Do not attempt to rinse and reinsert the device after initial insertion.
- Hold the surgical towel or 4 x 4 gauze pad away from the inflow and outflow windows, when controlling blood splatter during insertion of the Impella Heart Pump through the introducer.

If you miss a step in the process of setting up the Impella Catheter, or if you exceed the amount of time allowed to complete a step, the Automated Impella Controller will display additional instruction screens for continuing the setup process for some steps.

For the steps discussed below, All shaded steps require sterile technique.


Figure 1.4

### CASE START

1. Press the **START NEW CASE** soft button from the startup screen or plug in a new Impella Catheter. "Case Start" can also be selected by pressing the **MENU** soft key.
2. The controller displays the screen shown above.

### INSERT PURGE CASSETTE

1. Open the purge cassette package onto the sterile field.
2. If included, remove and discard the Y-connector and secure the yellow luer connector on the purge tubing to the sterile field.
3. Pass the purge cassette and spike off the sterile field.
4. Spike the glucose bag/bottle.
5. Press the **NEXT** soft button to continue.
6. Open the purge cassette door by pressing the release on the left side of the controller. Insert the purge cassette into the Automated Impella Controller.
7. Insert the purge cassette into the compartment on the front of the controller. Follow the diagram on the inside of the purge cassette door for proper placement.
8. Slide the purge disc into the slot to the right of the purge cassette until it snaps into place. The controller will automatically begin priming the purge cassette. Note: The instructions for inserting the purge disc appear if the disc is not snapped into place within 3 seconds of inserting the purge disc.

9. Extend the purge tubing and close the purge cassette door. There is sufficient room around the edges of the purge cassette door so that it will not pinch the purge tubing as it exits.
10. The controller automatically begins priming the purge cassette after it is inserted. The progress bar marks the progress of the purge cassette priming.

## CONNECT THE IMPELLA CATHETER

1. Remove the Impella Catheter from its package using sterile technique and inspect the catheter for damage.
2. Slide repositioning sheath towards the red pump handle.
3. Inspect the cable for damage, including damage to the connector pins at the controller end.
4. Pass the sterile connector cable from the Impella Catheter off the sterile field.


Figure 1.5

5. Open the cover on blue catheter plug by rotating clockwise. Line up the notch on the connector cable with the notch in the blue catheter plug on the front of the Automated Impella Controller and plug the cable into the controller.
6. Snap the purge clip (located on the pressure reservoir of the clear sidearm) to the connector cable as shown in Figure 1.5. **Note:** Snapping the plastic clip on the pressure reservoir to the connector cable is important to prevent the tube from kinking.
7. Once the purge cassette is primed and the controller detects that the connector cable is plugged in, it prompts you to connect the yellow luer to the Impella Catheter.


Figure 1.6

8. If included, remove and discard the Y-connector. Connect and tighten the yellow luer on the purge tubing to the Impella Catheter sidearm, if not already complete.


Figure 1.7

9. When the controller detects that the luer is connected, it automatically begins priming the purge lumen.
10. Once the purge cassette is primed and the luer is connected, the controller automatically advances to the next screen.
11. The first step on the next screen prompts you to enter the purge fluid information.

## ENTER PURGE FLUID DATA

1. Enter the purge fluid information.


Figure 1.8

2. To select the default values displayed on the screen, press the **ACCEPT** soft button. This will select those values and automatically advance to the next screen. **Note:** The default purge fluid values will be the purge fluid values from the last Case Start performed on the Automated Impella Controller.
3. To change the purge fluid information, press the **EDIT** soft buttons, scroll to the appropriate item and push the selector knob to select it or use the white arrow soft buttons. Then scroll through the values and push the selector knob or press **SELECT** to make a new selection. Press the **DONE** button to finish editing. The controller will use the default values if no other selections are made.
  - Purge fluid can be set to 50 mL, 100 mL, 250 mL, 500 mL, or 1000 mL.
  - Dextrose (Glucose) concentration can be set to 5%, 10% or 20%.
  - Heparin concentration can be set to 0 IU/ml, 5 IU/ml, 6.25 IU/ml, 10 IU/ml, 12.5 IU/ml, 20 IU/ml, 25 IU/ml, 40 IU/ml, 50 IU/ml.

## INSERTING THE IMPELLA CATHETER

**NOTE – Proper surgical procedures and techniques are the responsibility of the medical professional. The described procedure is furnished for information purposes only. Each physician must evaluate the appropriateness of the procedure based on his or her medical training and experience, the type of procedure, and the type of systems used.**

If the Impella Catheter is used in the OR as part of open heart surgery, manipulation may be performed only at the access site. Direct manipulation of the catheter assembly through the aorta or ventricle may result in serious damage to the Impella Catheter and serious injury to the patient.

When using a pigtail diagnostic catheter with side holes, ensure that the guidewire exits the end of the catheter and not the side hole. To do so, magnify the area one to two times as the guidewire begins to exit the pigtail.

During placement of the Impella Catheter, take care to avoid damage to the inlet area while holding the catheter and loading the placement guidewire.

1. Confirm purge fluid is exiting the Impella Catheter.
2. Obtain access to the femoral artery.
3. Insert a 5–8 Fr dilator over the 0.035 guidewire (provided) to predilate the vessel.
4. Remove the 5–8 Fr dilator over the 0.035 guidewire. Predilate the artery with the appropriate dilator prior to inserting the sheath. While inserting the introducer, hold the shaft of the introducer to slide it into the artery.

Figure 1.9

5. Administer heparin. When the ACT is greater than or equal to 250 seconds, remove the dilator.
6. Prior to catheter insertion, flush the introducer
7. Insert a diagnostic catheter (Abiomed recommends a 6 Fr AL1 or Multipurpose without side holes or 4–5 Fr pigtail with or without side holes) over a 0.035 inch diagnostic guidewire into the introducer and advance it into the left ventricle.

Figure 1.10

8. Remove the 0.035 inch diagnostic guidewire, leaving the diagnostic catheter in the ventricle. Form a curve or bend on the end of the 0.018 inch, 260 cm placement guidewire. Place the shaping tool just distal to the weld separating the shaping ribbon from the body of the placement guidewire. Bend the shaping ribbon against the tool, using minimal force. Do NOT use a shaping tool with a sharp tip or edge. Do NOT pull the shaping tool along the length of the shaping ribbon as this could strip the coil off the guidewire and cause it to unfurl and separate. Inspect the coil and guidewire for damage after shaping and before using.
9. Advance the placement guidewire into the apex of the left ventricle.
10. Remove the diagnostic catheter.
11. Slide the repositioning sheath towards the red pump plug prior to insertion.

### To backload the catheter using the EasyGuide lumen

12. Insert the placement guidewire into the red EasyGuide lumen at the tip of the pigtail. **(Once you remove the EasyGuide lumen from the Impella Catheter, do not attempt to reinsert it. If the catheter does not have a red EasyGuide lumen, follow the procedure outlined in step 13.)**
  - a. Advance the guidewire until it exits the red lumen near the label.
  - b. Remove the EasyGuide lumen by gently pulling the label while holding the Impella Catheter.
  - c. If you suspect that a portion of the red lumen remains in the catheter, do NOT use the Impella Catheter. Measure red lumen length using catheter markings (intact length is between 21.5 cm and 22.5 cm).
  - d. Skip to step 14 if the catheter is successfully backloaded on the guidewire.

Figure 1.11

### To backload the catheter without the EasyGuide lumen

13. Wet the cannula with sterile water and backload the catheter onto the placement guidewire. One or two people can load the catheter on the guidewire.

#### One-person technique

- a. Advance the guidewire into the Impella Catheter and stabilize the cannula between the fingers as shown in figure above. This prevents pinching of the inlet area. The guidewire must exit the outlet area on the inner radius of the cannula and align with the straight black line on the catheter as shown in figure below. The cannula can be hyperextended as necessary to ensure the guidewire exits on the inner radius of the cannula.

## Two-person technique

- b.** The scrub assistant can help stabilize the catheter by holding the catheter proximal to the motor. This will allow the implanting physician to visualize the inner radius. The guidewire must exit the outlet area on the inner radius of the cannula and align with the straight black line on the catheter, as shown in figure below. The physician can focus on advancing the guidewire and, if the cannula needs to be hyperextended, the scrub assistant is available to assist.


Figure 1.12

14. Advance the catheter through the hemostatic valve into the femoral artery and along the placement guidewire and across the aortic valve using a fixed-wire technique. Follow the catheter under fluoroscopy as it is advanced across the aortic valve, positioning the inlet area of the catheter 3.5 cm below the aortic valve annulus and in the middle of the ventricular chamber, free from the mitral valve chordae. Be careful not to coil the guidewire in the left ventricle.
15. Remove the placement guidewire.
16. Confirm position with fluoroscopy and confirm that an aortic waveform is displayed on the Automated Impella Controller.

## POSITIONING AND STARTING THE IMPELLA CATHETER

When the Impella Catheter is not correctly placed, there is no effective unloading of the ventricle. The patient may not be benefiting from the flow rate shown on the controller.

1. Reconfirm that the placement guidewire has been removed. Also reconfirm that the controller displays an aortic waveform and the radiopaque marker band is located at the aortic valve.
2. Press the **START IMPELLA** soft button. Press the **YES** soft button to confirm after wire is removed. The Impella will start in **AUTO** and automatically increase the flow rate over 30 seconds.
3. Once the controller has begun to run in **AUTO**, pressing the **FLOW CONTROL** soft button again opens the **FLOW CONTROL** menu with options for **AUTO**, and P-levels ranging from P-0 to P-9.
4. Wait 30 seconds for flow to reach its maximum value, then confirm correct and stable placement. Evaluate the catheter position in the aortic arch and remove any excess slack. The catheter should align against the lesser curvature of the aorta rather than the greater curvature. Verify placement with fluoroscopy and with the placement screen.
5. Reposition the catheter as necessary under fluoroscopy.

## MODES OF OPERATION

### AUTO

In **AUTO**, the Automated Impella Controller sets the motor speed of the Impella Catheter to achieve the maximum possible flow without causing suction. During **AUTO**, the highest the P-level the Impella CP with SmartAssist can run is P-9. After 3 hours of operation, the controller automatically switches to P-level mode. Upon transfer from **AUTO** mode to P-level mode, the controller displays the message shown below and the **AUTO** setting is no longer an option.


Figure 1.13 Transfer to P-level Mode

### P-LEVEL

In **P-LEVEL** mode you can select one of ten P-levels (P-0 to P-9) for the Impella Catheter (see Table below). Select the lowest P-level (P-2 or higher) that will enable you to achieve the flow rate necessary for patient support. At P-0, the Impella Catheter motor is stopped. If the Impella Catheter minimum flow is below 0.1 L/min then the controller will increase the motor speed to prevent retrograde flow.

Table 1.1 P-level Flow Rates for the Impella CP with SmartAssist Catheter

P-level	*Flow Rate (L/min)	Revolutions Per Minute (rpm)
P-0	0.0	0
P-1	0.0 - 0.9	23,000
P-2	1.1 - 2.1	31,000
P-3	1.6 - 2.3	33,000
P-4	2.0 - 2.5	35,000
P-5	2.3 - 2.7	37,000
P-6	2.5 - 2.9	39,000
P-7	2.9 - 3.3	42,000
P-8	3.1 - 3.4	44,000
P-9**	3.3 - 3.7	46,000

\*Flow rate can vary due to suction or incorrect positioning.

\*\* The peak flow rate at systole is up to 4.3 L/min at P-9

To operate the Impella Catheter in P-level mode:

1. Press the **FLOW CONTROL** soft button to open the **FLOW CONTROL** menu.
2. Turn the selector knob to increase or decrease the flow rate.
3. Press the selector knob to select the new flow rate.


Figure 1.14 Adjusting P-level

## ADJUSTING THE LV PLACEMENT SIGNAL

### CAUTIONS

- !** The LV placement signal and LV estimates are not displayed when pumps are running at P-3 or lower. Increase pump speed to P-4 or higher to re-enable signal
- !** Alarm conditions and low pump speeds may affect the LV placement signal and estimates
- !** Disruption of the outlet pressure placement signal, including alarms related to the Ao placement signal, will prevent calculation and display of an LV estimate. An operational Ao placement signal is required for the LV estimate
- !** LV placement signal calibration is not available if the P-level is less than P-4 or if Suction, Placement Signal Not Reliable, or positioning alarms are active
- !** Abnormal conditions, including cardiac arrhythmia, Ao-LV uncoupling, or aortic stenosis may limit the utility of the LV adjustment tool

### ADJUST LV PLACEMENT SIGNAL

Adjust the LV placement signal to reduce potential measurement variability. Adjust the LV Placement Signal when the white notification appears on the screen. The notification appears first when a suggested adjustment is calculated. If a suggested adjustment is calculated, then a second notification will appear after 24 hours of pump use.


Figure 1.15 White Notification for adjust LV

To adjust the LV Placement Signal:

1. Press the **MENU** soft button.
2. Select "Adjust LV Signal" option with rotary knob
3. Adjust the waveform using toggle arrows or rotary knob. The LV Adjustment defaults to the suggested adjustment value.
4. Press **DONE** to confirm suggested adjustment. LV waveform adjustments occur in increments of 1 mmHg from -60 to 60 mmHg. Recommended: Do not adjust to a value other than the suggested value.


Figure 1.16 Adjustment Tool


Figure 1.17 Placement Screen Post Adjustment

**Note:** It is not atypical to see a LV systolic value higher than the Ao systolic

## ENTERING CARDIAC OUTPUT


The cardiac output and cardiac power output metrics are for informational purposes only. Do not use for diagnostic purposes or patient monitoring. Independently verify all parameters displayed with a cleared or approved diagnostic device.

Enter Cardiac output into the Automated Impella Controller from a reference device, such as a Swan-Ganz catheter. Enter a new cardiac output every 8 hours. After 7 hours from entry, a white notification will trigger to enter a new cardiac output.

To enter the Cardiac Output:


1. Press **MENU** soft button.
2. Select "Enter Cardiac Output" using the rotary knob or the up/down arrow soft buttons.

3. Enter the total Cardiac Output, which can be any value from 0.0 to 10.0 L/min in increments of 0.1 L/min.
4. Press the **DONE** soft button or the rotary knob to complete.

**Figure 1.18 Enter Cardiac Output**

A Cardiac Output Confirmation will be displayed if a total Cardiac Output is entered that is less than or equal to the current Impella Flow (Figure 1.19). The **CONFIRM** soft button will use the current Impella Flow as the total Cardiac Output entry and the Native Cardiac Output will not be trended. The **BACK** soft button will go back to the Enter Cardiac Output screen. The **CANCEL** soft button will exit the workflow.

**Figure 1.19 Cardiac Output Confirmation**

#### CARDIAC OUTPUT, NATIVE CARDIAC OUTPUT & CARDIAC POWER OUTPUT CALCULATIONS

Once a cardiac output is entered, the Automated Impella Controller can calculate an initial cardiac power output and native cardiac output using the following equations:

$$\text{CPO} = (\text{CO} \times \text{MAP}) / 451$$

$$\text{NCO} = \text{CO} - \text{Impella Flow}$$

The Native Cardiac Output estimate is derived from a relationship between native function and aortic pulse pressure (PP). This relationship is linear and scaled by a calibration factor,  $\beta$  which may vary between patients and as an individual patient's condition changes. This relationships can be shown by the following equation:

$$\text{NCO} = \beta \times \text{PP}$$

Once the calibration factor is obtained, the Automated Impella Controller will continue to calculate the cardiac output, native cardiac output, and cardiac power output for the next 8 hours using the above equations.

If cardiac power output values are below or equal to 0.6, the value will display as yellow.

**Note:** Do not use values as a clinical diagnostic tool, this is for informational purposes only.

#### CARDIAC OUTPUT ENTRY REMINDERS

A white reminder to enter the CO will appear every 15 minutes for the last hour prior to the 8 hour timeout. If a new CO is not entered after 8 hours, the values will show as dash marks until a new cardiac output is entered.

A white reminder will also be displayed to update the CO if the Automated Impella Controller detects a significant change in the vascular state. This notification will be triggered if the average NCO or PP diverges from their initial values by a significant amount.

### USE OF THE REPOSITIONING SHEATH AND THE PEEL-AWAY INTRODUCER

1. Slide the repositioning sheath back to the red Impella plug.  
**Note:** For the 14 Fr x 25cm peel-away introducer, it may be necessary to pull the introducer over the repositioning sheath to remove the peel-away completely from the artery.
2. Remove the peel-away introducer completely from the artery over the catheter shaft to prevent trauma and significant bleeding and apply manual pressure above the puncture site.
3. Grasp the two "wings" and bend back until the valve assembly comes apart. Continue to peel the two wings until the introducer is completely separated from the catheter shaft.
4. Slide the repositioning sheath over the catheter shaft and advance it into the artery to the blue suture pads.
5. Secure the repositioning unit to the patient with the blue suture pads or a StatLock® stabilization device.
6. Evaluate the catheter position in the aortic arch and remove any excess slack. The catheter should align against the lesser curvature of the aorta rather than the greater curvature. Verify placement with fluoroscopy and with the placement signal.
7. Attach the anticontamination sleeve to the blue section of the repositioning sheath. Lock the anchoring ring in place by turning it clockwise. Secure the catheter shaft in place by tightening the connected anchoring ring.
8. Carefully extend the anticontamination sleeve to maximum length and secure the end closest to the red Impella plug by tightening the anchoring ring.

### PATIENT WEANING

Weaning the patient from the Impella Catheter is at the discretion of the physician. The LVEDP/CO Trend can be used while weaning. Refer to the Automated Impella Controller Instructions for Use document for more information.

***The following weaning instructions are provided as guidance only.***

1. To initiate weaning, press **FLOW CONTROL** and reduce P-level by 2 level increments over time intervals as cardiac function allows.
2. Keep Impella Catheter P-level at P-2 or above until the catheter is ready to be explanted from the left ventricle.
3. When the patient's hemodynamics are stable, reduce the P-level to P-2 and pull the Impella Catheter back across the aortic valve into the aorta.
4. If the patient's hemodynamics remain stable, follow instructions in the next section for removing the Impella Catheter.

## REMOVING THE IMPELLA CATHETER

The Impella Catheter can be removed after weaning when the introducer is still in place or when the catheter is secured with the repositioning sheath.

Removal of the Impella Catheter must be completed with care to avoid damage to the catheter assembly. Do NOT attempt to remove the Impella Catheter through the repositioning sheath. The Impella Catheter will not pass through the repositioning sheath.

### REMOVING THE IMPELLA CATHETER WITH THE INTRODUCER IN PLACE

1. Reduce the P-level to P-0.
2. Remove the Impella Catheter through the introducer.
3. Wait until ACT drops below 150 seconds.
4. When ACT is below 150 seconds, remove the introducer.
5. Disconnect the connector cable from the Automated Impella Controller and turn the controller off by pressing the power switch on the side of the controller for 3 seconds.
6. Apply manual compression for 40 minutes or per hospital protocol.

### REMOVING THE IMPELLA CATHETER SECURED WITH THE REPOSITIONING SHEATH

1. When ACT is below 150 seconds, press **FLOW CONTROL** and reduce the P-level to P-0.
2. Remove the Impella Catheter and repositioning sheath together (the catheter will not come through the repositioning sheath).
3. Disconnect the connector cable from the Automated Impella Controller and turn the controller off by pressing the power switch on the side of the controller for 3 seconds.
4. Apply manual compression for 40 minutes or per hospital protocol.

### REMOVING THE IMPELLA CP WITH SMARTASSIST CATHETER WHILE MAINTAINING GUIDEWIRE ACCESS

1. Remove the stylet from the guidewire access port.


Figure 1.20 Removing the Stylet

2. Aspirate using a syringe to ensure that the line is patent; confirm pulsatile blood flow from entrance port.
3. Advance the 0.035" (or smaller) guidewire with an atraumatic tip through the guidewire access port using the supplied cheater.


Figure 1.21 Inserting the Guidewire with the Cheater

4. Advance the guidewire tip into the descending aorta under fluoroscopic guidance.
5. Loosen the anticontamination sleeve by rotating the Tuohy-Borst valve counterclockwise.


Figure 1.22 Loosening the Tuohy-Borst Valve

6. Anchor the guidewire and repositioning unit and withdraw the catheter until the distal end of the catheter reaches the distal tip of the repositioning unit.
7. With the guidewire anchored and pressure applied to the access site, completely remove the Impella Catheter and repositioning sheath together.

## OPERATING THE IMPELLA CATHETER WITHOUT HEPARIN IN THE PURGE SOLUTION

The Impella Catheter is designed to be operated with a purge solution that contains heparin. Operation of the system without heparin in the purge solution has not been tested. In the event that a patient is intolerant to heparin, due to heparin-induced thrombocytopenia (HIT) or bleeding, physicians should use their clinical judgment to assess the risks versus benefits of operating the Impella System without heparin.

If it is in the best interest of the patient to operate the system without heparin, the glucose solution is still required, and physicians should consider *systemic delivery* of an alternative anticoagulant. DO NOT add any alternative anticoagulant (such as a direct thrombin inhibitor) to the purge fluid. The Impella Catheter has not been tested with any alternative anticoagulants in the purge solution.

## ALARMS OVERVIEW

For a full list of all Impella alarms, refer to the Automated Impella Controller Instructions for Use document.

The Automated Impella Controller monitors various functions to determine whether specific operational parameters are within expected limits. When a parameter goes outside of its specified limits, the controller sounds an alarm tone and displays an alarm message that can be viewed on the display screen on the front of the controller. The alarm tone indicates the severity of the alarm. The alarm message on the display screen is color-coded for severity and provides details on the cause of the alarm and how to resolve the alarm. After muting an alarm, if another alarm occurs it will only be heard and displayed if it is a higher priority alarm than the one that was muted.

Category	Description	Audible Indicator*	Visual Indicator
Advisory	Notification	1 beep every 5 minutes	Alarm header on white background
Serious	Abnormal situation. Prompt action needed	3 beeps every 15 seconds	Alarm header on yellow background
Critical	High priority. Immediate action needed.	10 beeps every 6.7 seconds	Alarm header on red background

\* Sound pressure of audible alarm indicators is 60-80 dBA at a measurement distance of 1 m.

Table 1.2 Alarm Levels

### MUTE ALARM FUNCTION

Pressing the **MUTE ALARM** button on the upper right of the Automated Impella Controller display screen will silence the audible alarm indicator for 2 minutes (for red or yellow alarms) or 5 minutes (for white advisory alarms). When an alarm is silenced, the words “**MUTE ALARM**” next to the button are replaced by the mute alarm indicator, a crossed-out bell icon.

The audible indicator will shut off if an alarm condition is resolved before you press **MUTE ALARM**. The visual message, however, will continue to be displayed, with the alarm header on a gray background, for 20 minutes or until you press **MUTE ALARM**. This allows you to identify the alarm that occurred.

## SYMBOLS

Table 1.3 Symbols

		Caution; follow instructions for use
		Defibrillator-proof type CF equipment
		Keep dry
		Storage temperature (eg, 10°C to 30°C)
		Declares conformity with Directive 93/42/EEC for medical devices or with Regulation (EU) 2017/745 on Medical Devices, and Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment
		Date of manufacture (eg, 2019-11-30)
		Protect from sunlight
		Symbol for lot designation; the manufacturer's lot designation must be stated after the LOT symbol
		Abiomed part number (eg, part number 123456)
		Manufacturer's serial number (eg, serial number 123456)
		Use-by date (eg, use before 2019-11-30)
		Do not reuse
		Sterilized using ethylene oxide
		Medical Device
		Federal (USA) law restricts this device to sale by or on order of a physician: Prescription only.
		Do not use if package is damaged
		Single sterile barrier system with protective packaging inside
	<b>Do Not Flush</b>	Do NOT flush with syringe
	<b>Glucose</b>	Use glucose in the purge fluid
<b>Do Not Clean with Alcohol</b>		Do not use alcohol or alcohol-based products for cleaning.

## IMPELLA CATHETER PARAMETERS

Table 1.4 Impella Catheter Parameters

Impella CP® with SmartAssist®	
<b>Speed range</b>	0 to 46,000 rpm
<b>Power consumption</b>	24 W
<b>Voltage</b>	Max. 20 V DC
<b>Maximum Mean Flow</b>	3.7 L/min
<b>Purging the Impella Catheter</b>	
Recommended purge fluid	5% glucose solution with heparin concentration of 50 IU per mL
Glucose concentration	5% to 20%
Purge pressure	300 to 1100 mmHg
Infusion rate	2 to 30 mL/h
<b>Maximum duration of use</b>	Up to 5 days
<b>Catheter dimensions</b>	
Length of invasive portion (without catheter)	150 ± 3 mm
Diameter	Max. 4.9 mm (nom. 4.7 mm)
<b>Cable length</b>	300 cm
<b>Classification per DIN EN 60601-1</b>	Protection class I, degree of protection: CF (Automated Impella Controller and Impella Catheter)
<b>Classification per Regulation (EU) 2017/745</b>	Class III
<b>Latex</b>	Not made with natural rubber latex

## PUMP METRICS SPECIFICATIONS

Table 1.5 Pump Metrics Specifications

Frequency	Range	Accuracy*
Pump Outlet (Aortic) Pressure	0 - 200 mmHg	4.2 mmHg
Pump Inlet (Left Ventricular) Pressure	0 to 200 mmHg	5.3 mmHg
Pulse Pressure	1 to 100 mmHg	3.5 mmHg
Native Cardiac Output	0 to 5.0 L/min.	0.4 L/min.
Cardiac Power Output	0 to 3.0 Watts	0.2 Watts

\* the measured root mean square error (of multiple measurements)

Metrics displayed are not intended for diagnostic use. All parameters displayed must be verified independently using either a cleared or approved diagnostic device, and must not be used for patient monitoring.

## IMPELLA CP® WITH SMARTASSIST® CATHETER DIMENSIONS


Figure 1.23 Impella CP with SmartAssist Catheter Dimensions

## DISPOSING OF THE IMPELLA CATHETER AND ACCESSORIES

The Impella Catheter is a disposable item that must be disposed of in accordance with hospital regulations for blood contaminated materials.

The Automated Impella Controller is marked according to Directive 2012/19/EU. Devices sold within the EEC can be returned to Abiomed Europe GmbH for correct disposal.

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Recovering hearts. Saving lives.

**Clinical support 24 hours per day, 7 days a week:**

**1-800-422-8666 (US)**  
**+00800 0 22 466 33 (EU)**  
**+0800 100 490 (France)**

**[www.abiomed.com](http://www.abiomed.com)**

**Abiomed, Inc.**  
22 Cherry Hill Drive  
Danvers, Massachusetts 01923 USA  
Voice: +1 978-646-1400  
Facsimile: +1 978-777-8411  
Email: clinical@abiomed.com

**Abiomed Europe GmbH**  
Neuenhofer Weg 3  
52074 Aachen, Germany  
Voice: +49 (0) 241 8860-0  
Facsimile: +49 (0) 241 8860-111  
Email: europe@abiomed.com

**CE 0482**

©2025 ABIOMED, Inc. All rights reserved.  
The ABIOMED logo and ABIOMED, Impella, Impella 2.5, Impella CP, Impella CP with SmartAssist, Impella 5.0, Impella 5.5 with SmartAssist, Impella LD, Impella RP, and Impella Connect are registered trademarks of Abiomed, Inc. in the U.S.A. and certain foreign countries.  
Impella® is a registered trademark of Abiomed Europe GmbH, a wholly owned subsidiary of Abiomed, Inc., in the U.S.A. and certain foreign countries.